

IAPSO MIO

Dominique Lefèvre (CNRS), dominique.lefevre@mio.osupytheas.fr
Anthony Bosse (Aix-Marseille Université), anthony.bosse@mio.osupytheas.fr

CTD and oxygen adjustment procedure

The general qualification of the CTD data follows the following workflow, where the final .cnv and .btl files contain adjusted salinity and dissolved oxygen data relative to the reference autosal and winkler data.




Figure 1: General workflow for the CTDO2 adjustments in delayed mode.

Seabird Processing

Folder architecture: a working directory is setup with the following folders:

SBE-profile-data/raw (containing raw data *.hex, *.bl, *.hdr)
SBE-profile-data/xmlcon (containing original config files *.xmlcon)
SBE-profile-data/cnv (where processed *.cnv files will be stored)
SBE-profile-data/btl (where *.btl will be stored)

SBE_Pprocess/ (with a batch executive file will check for the presence of a cnv file in /cnv, will run "sbebatch" if not, until all files are processed)

SBE_Process/psa_doubleCT (containing the *.psa configuration files for the data processing including traitementctd.txt which gives the order of steps to follow)

Example of the different steps done by SBE processing:

```
datcnv /i..\SBE-profile-data\raw\%1.hex /c..\SBE-profile-data\xmlcon_adj\%1.XMLCON /o..\SBE-profile-data\cnv_adj
/ppsa_doubleCT\datcnv.psa /f%1
wildedit /i..\SBE-profile-data\cnv_adj\%1.cnv /o..\SBE-profile-data\cnv_adj/ppsa_doubleCT\wildedit.psa /f%1
wildedit /i..\SBE-profile-data\cnv_adj\%1.cnv /o..\SBE-profile-data\cnv_adj/ppsa_doubleCT\wildedit.psa /f%1
filter /i..\SBE-profile-data\cnv_adj\%1.cnv /o..\SBE-profile-data\cnv_adj/ppsa_doubleCT\filter.psa /f%1
alignctd /i..\SBE-profile-data\cnv_adj\%1.cnv /o..\SBE-profile-data\cnv_adj/ppsa_doubleCT\alignctd.psa /f%1
celltm /i..\SBE-profile-data\cnv_adj\%1.cnv /o..\SBE-profile-data\cnv_adj/ppsa_doubleCT\cellctm.psa /f%1
loopedit /i..\SBE-profile-data\cnv_adj\%1.cnv /o..\SBE-profile-data\cnv_adj/ppsa_doubleCT\loopedit.psa /f%1
derive /i..\SBE-profile-data\cnv_adj\%1.cnv /c..\SBE-profile-data\xmlcon_adj\%1.XMLCON /o..\SBE-profile-data\cnv_adj
/ppsa_doubleCT\derive.psa /f%1
binavg /i..\SBE-profile-data\cnv_adj\%1.cnv /o..\SBE-profile-data\cnv_adj/ppsa_doubleCT\binavg.psa /f%1
bottlesum /i..\SBE-profile-data\cnv_adj\%1.ros /c..\SBE-profile-data\xmlcon_adj\%1.XMLCON /o..\SBE-profile-data\bt\adj
/ppsa_doubleCT\bottlesum.psa /f%1
```

The following configuration for each psa is used:

datcnv.psa

```
# datcnv_skipover = 0
# datcnv_ox_hysteresis_correction = yes
```

wildedit.psa

```
# wildedit_pass1_nstd = 2.0
# wildedit_pass2_nstd = 10.0
# wildedit_pass2_mindelta = 0.000e+000
# wildedit_npoint = 100
# wildedit_vars = prDM t090C t190C c0S/m c1S/m
# wildedit_excl_bad_scans = yes
```

filter.psa

```
# filter_low_pass_tc_A = 0.030
# filter_low_pass_tc_B = 0.150
# filter_low_pass_A_vars = turbWETntu0 spar wetStar
# filter_low_pass_B_vars = prDM
```

alignctd.psa

```
# alignctd_adv = sbeox0V 2.000
```

celltm.psa

```
# celltm_alpha = 0.0300, 0.0300
# celltm_tau = 7.0000, 7.0000
# celltm_temp_sensor_use_for_cond = primary, secondary
```

loopedit.psa

```
# loopedit_minVelocity = 0.100
# loopedit_surfaceSoak: minDepth = 6.0, maxDepth = 25, useDeckPress = 1
# loopedit_excl_bad_scans = yes
```

derive.psa

```
# derive_time_window_docdt = seconds: 2
# derive_ox_tau_correction = yes
```

binavg.psa

```
# binavg_bintype = decibars
# binavg_binsize = 1
# binavg_excl_bad_scans = yes
# binavg_skipover = 0
# binavg_omit = 0
```

To ease the work, “*.bti” files are merged into one single file using a R script, from which adjustments can be made.

Cf annexe

1. Raw data batch processing

As shown by the general workflow, the first step consists in a batch processing of the raw CTD files, producing .cnv and .blt files.

2. Conductivity-Salinity validation using Salinity bottles

The conductivity correction is deduced from the following steps (see also AUTOSAL_comparison.ipynb).

First, Autosal salinities are plotted against bottle salinities and the presence of outliers stands out clearly (Figure 1).






Figure 2: (left) Autosal and bottle salinities versus pressure. (right) Difference of Autosal and bottle salinities against pressure for the two CTD sensors.

A subsample of the salinities to be used is determined by using a standard Interquartile range (IQR) method excluding data outside the interval $[Q1 - 1.5 \cdot (Q3 - Q1) ; Q3 + 1.5 \cdot (Q3 - Q1)]$ with Q1 and Q3 the 25th and 75th percentile of the data distribution. This reduces the

number of samples from 1153 to 1003 and 986, for CTD1 and CTD2 respectively (Figure 3-left).




The distribution of differences can be plotted and examined. It shows small medians of -0.003 and -0.006 psu respectively with CTD1 and CTD2 (Figure 4).



Note that there seems to be a shift in the distribution with pressure, especially for CTD1. The order of magnitude of this shift is not larger than about 0.001 psu per 1000 dbar, which

remains small compared to the variability of the signal. The correction applied here does not correct for any pressure dependence. This might have to be corrected in a second step in order to get more accurate salinity data below 2000 dbar.



Salinity differences are also plotted against time, to check that no temporal drift was present. There is a general decrease in salinity difference of about -0.003psu in the 3 weeks of the sampling.

The next step is the determination of the slope coefficient. For a set of samples, the slope coefficient is calculated as

(https://www.seabird.com/cms-portals/seabird_com/cms/documents/training/Module10_DataAccuracyFieldCals.pdf):

$$\text{slope} = \frac{\sum_{i=1}^n \alpha_i \beta_i}{\sum_{i=1}^n \alpha_i \alpha_i}$$

where : n = number of samples, alpha = CTD conductivity, beta = true (bottle sample) conductivity

A general coefficient is calculated for each CTD. In our case, it gives:

Slope correction for CTD1 | CTD2 : 0.999995 | 0.999986

The first option is to set those slope coefficients constant for CTD1 and CTD2 in order to correct the whole data set. But, as shown by Figure 5, the salinity difference seems to vary with time. So we decided to retain a time dependence in the slope coefficients (Figure 6). This means that for each CTD cast, a slope correction coefficient will be interpolated from the slope/time relationship defined by the 5-day rolling median signal of Figure 6.




Figure 6: Slope coefficient for the two CTDs according to time. The scattered dots represent slope coefficient computed for each single cast. A 5-day rolling median is then applied to deduce a time-varying slope correction coefficient.

Before reprocessing the whole dataset, conductivity is recalculated using the inferred slope coefficient. The results show better with a median value centered (0.00009/0.00008psu, with standard deviation of 0.0019psu) and no more time variation in salinity difference as expected. (Figure 7 and 8)

Note that the initial difference in salinity was already before correction within the range of precision of the CTD instrument (0.0001 +/- 0.002 psu) and the present correction does not massively increase the precision of the final data set in terms of absolute salinity values.




Figure 5: Autosal versus bottle salinity difference as a function of time after correction with a time-varying slope coefficient. The dashed lines show the rolling median over a 5-day window.




Figure 4: Histogram of the difference between Autosal and bottle salinities after correction with a time-varying slope coefficient.

3. Oxygen SBE43 configuration adjustment using discrete Winkler samples.

From the merged file (BTL Winkler and SBE43), use the Seabird application note 64.2 (http://ftp.seabird.com/application_notes/AN64-2.htm) to adjust SOC, Offset and E parameter in the following equation:

Based on Seabird algorithm

$$OX = Soc * [V + Voffset] * OxSol(Theta, S) * (1.0 + A * T90C + B * T90C^2 + C * T90C^3) * \exp(E * P / K)$$

Ox, Oxygen concentration in $\mu\text{mol kg}^{-1}$

P : Pressure in db

T90, in situ temperature $^{\circ}\text{C}$

Theta, potential temperature $^{\circ}\text{C}$

S, salinity

V, SBE43 O2 voltage in volt

$OxSol(Theta, Salinity) = \text{Garcia_Gordon1992_Benson_Krause coefficients } \mu\text{mol kg}^{-1}$

*It is to be noticed that there is a mixed usage of in situ and potential temperature in this equation !

Adjust SBE_Oxygen value (O2_SBE43) is calculated based on this equation for each sample where there is a O2_Winkler sample (O2_W)

The sum of the square distance between each O2_W and O2_SBE_43 is computed.

This sum is then minimised using the evolutionary Excel ® solver function by adjusting the SOC, Offset and E parameters of the Seabird ® equation

For each sample, the value of the distance between the O2W and O2SBE43 is computed and any distance greater than $\pm 2 * Sd$ is removed.


The solver function is run again until there are no more distance value greater than $\pm 2 * Sd$.

It can take several iterations, in this case it took 10 iterations.


The adjusted SOC, Offset and E value are then used in the xlmcon file to generate the new adjusted profiles.

To do this we use Excel® spreadsheet and VBA code under Excel provided in the spreadsheet.


This process is summarised on the following flow chart




Using this iterative process, 212 points have been removed. 921 samples were kept for the adjustment of the SBE43 with winkler samples.




These 3 charts are showing the distribution of the removed data points (Winkler samples) to assess any potential bias such as faulty Niskin bottles. These points are randomly distributed and do not exhibit any obvious bias.





Plot property of the SBE43 vs Winkler samples before adjustment and after adjustment. The effect of “E” pressure coefficient can be identified from the disappearance of the scattering.



Distribution of the residuals before and after adjustment against Niskin bottle. No Bias can be observed.



Distribution of the residuals before and after adjustment against cast number.
No time drift can be derived from this plot.



Vertical distribution of O2 concentration (left) vs Pressure of the BTL-SBE43 and Winkler samples and O2 residuals (right). residuals are homogeneously distributed around the zero

value throughout the water column, providing some confidence in the adjustment procedure of the coefficients.

The adjusted coefficient are

	Original	Adjusted
SOC	5.5183E-01	5.7136E-01
Offset	-5.0930E-01	-5.0233E-01
E	3.6000E-02	3.7649E-02

RMSE = 1.28

4. Reprocessing of CTD files with adjusted configuration files

A final step makes sure that the correction have well been applied using the adjusted *_adj.xmlcon files.




Figure : Overview of the reprocessed CTD files (temperature, salinity and oxygen) in the top 200m.




Figure : Overview of the reprocessed CTD files (temperature, salinity and oxygen) with full-depth y-axis.

Annexe Code R

Code R to merge single *.btl file into one file to ease the adjustment process.

```
##IAPSO
####HOMEPATH../IAPSO/SBE-profile-data/btl/" ### To be ADJUSTED for each user ####

file_label<-list.files(HOMEPATH)
source("../IAPSO/CODE/read_btl_info_V_IAPSO.R") # call function
source("../IAPSO/CODE/read_btl_V_IAPSO.R") # Call function
#      Bottle      Date Potemp090C      Sal00 Potemp190C      Sal11 Sbeox0ML/L
Sbox0Mm/Kg OxsolMm/Kg      PrDM      T090C      T190C C0S/m      C1S/m      Sbeox0V
      Scan
#header<-c("Bottle","Month","Day","Year","Latitude","Longitude","PrDM","T090C","C0S/m","Sal00","T1
90C","C1S/m","Sal11","FIC","CStarTr0","Sbeox0V","Sbox0Mm/Kg")

header<-c("Bottle","Month","Day","Year",      "Potemp090C",      "Sal00", "Potemp190C",
"Sal11", "Sbeox0ML/L", "Sbox0Mm/Kg", "OxsolMm/Kg",      "PrDM",      "T090C",      "T190C",
      "C0S/m" , "C1S/m" , "Sbeox0V",      "Scan")

header_Out<-c("Bottle","Month","Day","Year","Latitude","Longitude",      "Potemp090C",
"Sal00", "Potemp190C",      "Sal11", "Sbeox0ML/L", "Sbox0Mm/Kg", "OxsolMm/Kg",      "PrDM",
"T090C",      "T190C",      "C0S/m" ,      "C1S/m" , "Sbeox0V","BTLFILE")
#####
for(j in 1:length(file_label)){
#####
  dfbtl<-read_btl(btl_ori_file=paste0(HOMEPATH,file_label[j]),btl_header=header,btl_file_name =
file_label[j])
  if (j==1) dfbtl_Out<-dfbtl else dfbtl_Out<-rbind(dfbtl_Out,dfbtl)
}
dfbtl_Out<-dfbtl_Out[,header_Out]
write.csv(dfbtl_Out,file=" ../IAPSO/SBE-profile-data/IAPSO_BTL_adj.csv", row.names = FALSE)

##### FUNCTION READ_BTL
read_btl_info<-function(btl_ori_file,header_max=500){
#####
#####
line2head=0
line2skip2all=0
for (i in 1:header_max){

line<-scan(btl_ori_file,skip=i,nline=1,what="character",sep="",quiet=TRUE,nmax=1)
      if (length(line)==0) break
      if (line=="") line2head=i
      if (line=="Bottle") line2skip2all=i
    }
    Allline<-i
    n_btl<-(Allline-line2skip2all-2)/2
    print(paste("number of nskin:",n_btl))
  for (i in 1:line2head){

line<-scan(btl_ori_file,skip=i,nline=1,what="character",sep="",quiet=TRUE,nmax=3)
```

```

        if (line[2]=="NMEA" & line[3]=="Latitude") line2LAT<-i
        if (line[2]=="NMEA" & line[3]=="Longitude") line2LON<-i
        if (line[2]=="NMEA" & line[3]=="UTC") line2TIME<-i
    }
    if (exists("line2LAT")==FALSE & exists("line2LON")==FALSE & exists("line2TIME")==FALSE) print("no
    NMEA lat,lon,time available") else {
    lat<-scan(btl_ori_file,skip=line2LAT,nline=1,what="character",sep="=",quiet=TRUE,nmax=4)[2]
    lon<-scan(btl_ori_file,skip=line2LON,nline=1,what="character",sep="=",quiet=TRUE,nmax=4)[2]
    time<-scan(btl_ori_file,skip=line2TIME,nline=1,what="character",sep="=",quiet=TRUE,nmax=4)[2]
    print(paste("Latitude:",lat))
    print(paste("Longitude:",lon))
    print(paste("Time:",time))
    }
    #####
    #btl<-as.data.frame(matrix(NA,ncol=length(btl_header)))
    #####
    line2read<-seq(2,n_btl*2,2)
    line2skip<-line2skip2all+line2read
    headerinfo<-scan(file=btl_ori_file,skip=line2skip2all,nline=1,what="character",sep="",quiet=TRUE)
    print("header line is :")
    print(headerinfo)
    for (k in 1:length(line2read)) {
        line<-read.table(file=btl_ori_file,skip=line2skip[k],nrows=1)
        if (k==1) {
            print("First data line is :")
            print(line)
        }
    }
    #        btl[k,]<-line
    }

##### FUNCTION READ_BTL
read_btl_info<-function(btl_ori_file,header_max=500){
#####
#####
line2head=0
line2skip2all=0
for (i in 1:header_max){

line<-scan(btl_ori_file,skip=i,nline=1,what="character",sep="",quiet=TRUE,nmax=1)
        if (length(line)==0) break
        if (line=="") line2head=i
        if (line=="Bottle") line2skip2all=i
        }
        Allline<-i
        n_btl<-(Allline-line2skip2all-2)/2
        print(paste("number of nskin:",n_btl))
    for (i in 1:line2head){

line<-scan(btl_ori_file,skip=i,nline=1,what="character",sep="",quiet=TRUE,nmax=3)
        if (line[2]=="NMEA" & line[3]=="Latitude") line2LAT<-i
        if (line[2]=="NMEA" & line[3]=="Longitude") line2LON<-i

```



```

        if (line[2]=="NMEA" & line[3]=="UTC") line2TIME<-i
    }
if (exists("line2LAT")==FALSE & exists("line2LON")==FALSE & exists("line2TIME")==FALSE) print("no
NMEA lat,lon,time available") else {
lat<-scan(btl_ori_file,skip=line2LAT,nline=1,what="character",sep="=",quiet=TRUE,nmax=4)[2]
lon<-scan(btl_ori_file,skip=line2LON,nline=1,what="character",sep="=",quiet=TRUE,nmax=4)[2]
time<-scan(btl_ori_file,skip=line2TIME,nline=1,what="character",sep="=",quiet=TRUE,nmax=4)[2]
print(paste("Latitude:",lat))
print(paste("Longitude:",lon))
print(paste("Time:",time))
}
#####
#btl<-as.data.frame(matrix(NA,ncol=length(btl_header)))
#####
line2read<-seq(2,n_btl*2,2)
line2skip<-line2skip2all+line2read
headerinfo<-scan(file=btl_ori_file,skip=line2skip2all,nline=1,what="character",sep="",quiet=TRUE)
print("header line is :")
print(headerinfo)
for (k in 1:length(line2read)) {
    line<-read.table(file=btl_ori_file,skip=line2skip[k],nrows=1)
    if (k==1) {
        print("First data line is :")
        print(line)
    }
#        btl[k,]<-line
    }
}

```

#####

Code VBA Excel

#####

Sub MergeBTL_Chemistry_IAPSO()

```

File = "IAPSO_BTL_Merged.xlsm"
O2W_S = "CE17007BottleODV"
BTL = "IAPSO_BTL"
O2 = "Oxygen"
Sal = "Salinity"

```

i = 4 ' Index of F_O2

j = 3 ' Index of F_OUT

k = 2

M = 2

'Reading source of chemistry

For i = 2 To 1196

Cast = Workbooks(File).Worksheets(O2W_S).Cells(i, 1)

' Ensure Niskin are ranked othewise, test check NiskinNumber?

```

Niskin = Workbooks(File).Worksheets(O2W_S).Cells(i, 8)

' Check if Sal or O2W exist
O2W = Workbooks(File).Worksheets(O2W_S).Cells(i, 14)
Salinity = Workbooks(File).Worksheets(O2W_S).Cells(i, 12)
If O2W <> "" And O2W <> -999 Then
'get BTL info
For j = 2 To 1208
    CastBtl = Workbooks(File).Worksheets(BTL).Cells(j, 20)
    NiskinBtl = Workbooks(File).Worksheets(BTL).Cells(j, 1)
    If CastBtl = Cast And NiskinBtl = Niskin Then
        Workbooks(File).Worksheets(O2).Cells(k, 1) = Cast
        Workbooks(File).Worksheets(O2).Cells(k, 2) = Niskin
        Workbooks(File).Worksheets(O2).Cells(k, 3) = O2W
        'Get BTL values to complete
        '
        Temp = Theta
        '
        Salinity = Salinity
        '
        SolO2_Theta_molkg = SolO2μmolkg(Temp, Salinity)
        Workbooks(File).Worksheets(O2).Cells(k, 4) =
Workbooks(File).Worksheets(BTL).Cells(j, 2) 'Month
        Workbooks(File).Worksheets(O2).Cells(k, 5) =
Workbooks(File).Worksheets(BTL).Cells(j, 3) 'Day
        Workbooks(File).Worksheets(O2).Cells(k, 6) =
Workbooks(File).Worksheets(BTL).Cells(j, 4) 'Year
        Workbooks(File).Worksheets(O2).Cells(k, 7) =
Workbooks(File).Worksheets(BTL).Cells(j, 5) 'Latitude
        Workbooks(File).Worksheets(O2).Cells(k, 8) =
Workbooks(File).Worksheets(BTL).Cells(j, 6) 'Longitude
        Workbooks(File).Worksheets(O2).Cells(k, 9) =
Workbooks(File).Worksheets(BTL).Cells(j, 7) 'Potemp090C
        Workbooks(File).Worksheets(O2).Cells(k, 10) =
Workbooks(File).Worksheets(BTL).Cells(j, 8) 'Sal00
        Workbooks(File).Worksheets(O2).Cells(k, 11) =
Workbooks(File).Worksheets(BTL).Cells(j, 9) 'Potemp190C
        Workbooks(File).Worksheets(O2).Cells(k, 12) =
Workbooks(File).Worksheets(BTL).Cells(j, 10) 'Sal11
        Workbooks(File).Worksheets(O2).Cells(k, 13) =
Workbooks(File).Worksheets(BTL).Cells(j, 11) 'Sbeox0ML/L
        Workbooks(File).Worksheets(O2).Cells(k, 14) =
Workbooks(File).Worksheets(BTL).Cells(j, 12) 'Sbox0Mm/Kg
        Workbooks(File).Worksheets(O2).Cells(k, 15) =
Workbooks(File).Worksheets(BTL).Cells(j, 13) 'OxsolMm/Kg
        Workbooks(File).Worksheets(O2).Cells(k, 16) =
Workbooks(File).Worksheets(BTL).Cells(j, 14) 'PrDM
        Workbooks(File).Worksheets(O2).Cells(k, 17) =
Workbooks(File).Worksheets(BTL).Cells(j, 15) 'T090C
        Workbooks(File).Worksheets(O2).Cells(k, 18) =
Workbooks(File).Worksheets(BTL).Cells(j, 16) 'T190C

```

```

        Workbooks(File).Worksheets(O2).Cells(k, 19) =
Workbooks(File).Worksheets(BTL).Cells(j, 17) 'C0S/m
        Workbooks(File).Worksheets(O2).Cells(k, 20) =
Workbooks(File).Worksheets(BTL).Cells(j, 18) 'C1S/m
        Workbooks(File).Worksheets(O2).Cells(k, 21) =
Workbooks(File).Worksheets(BTL).Cells(j, 19) 'Sbeox0V
        k = k + 1
    End If

```

```

Next j
End If
If Salinity <> "" And Salinity <> -999 Then
'get BTL info
For j = 2 To 1208
    CastBtl = Workbooks(File).Worksheets(BTL).Cells(j, 20)
    NiskinBtl = Workbooks(File).Worksheets(BTL).Cells(j, 1)
    If CastBtl = Cast And NiskinBtl = Niskin Then
        Workbooks(File).Worksheets(Sal).Cells(M, 1) = Cast
        Workbooks(File).Worksheets(Sal).Cells(M, 2) = Niskin
        Workbooks(File).Worksheets(Sal).Cells(M, 3) = Salinity
        'Get BTL values to complete
        '
        Temp = Theta
        '
        Salinity = Salinity
        '
        SolO2_Theta_molkg = SolO2μmolkg(Temp, Salinity)
        Workbooks(File).Worksheets(Sal).Cells(M, 4) =
Workbooks(File).Worksheets(BTL).Cells(j, 2) 'Month
        Workbooks(File).Worksheets(Sal).Cells(M, 5) =
Workbooks(File).Worksheets(BTL).Cells(j, 3) 'Day
        Workbooks(File).Worksheets(Sal).Cells(M, 6) =
Workbooks(File).Worksheets(BTL).Cells(j, 4) 'Year
        Workbooks(File).Worksheets(Sal).Cells(M, 7) =
Workbooks(File).Worksheets(BTL).Cells(j, 5) 'Latitude
        Workbooks(File).Worksheets(Sal).Cells(M, 8) =
Workbooks(File).Worksheets(BTL).Cells(j, 6) 'Longitude
        Workbooks(File).Worksheets(Sal).Cells(M, 9) =
Workbooks(File).Worksheets(BTL).Cells(j, 7) 'Potemp090C
        Workbooks(File).Worksheets(Sal).Cells(M, 10) =
Workbooks(File).Worksheets(BTL).Cells(j, 8) 'Sal00
        Workbooks(File).Worksheets(Sal).Cells(M, 11) =
Workbooks(File).Worksheets(BTL).Cells(j, 9) 'Potemp190C
        Workbooks(File).Worksheets(Sal).Cells(M, 12) =
Workbooks(File).Worksheets(BTL).Cells(j, 10) 'Sal11
        Workbooks(File).Worksheets(Sal).Cells(M, 13) =
Workbooks(File).Worksheets(BTL).Cells(j, 11) 'Sbeox0ML/L
        Workbooks(File).Worksheets(Sal).Cells(M, 14) =
Workbooks(File).Worksheets(BTL).Cells(j, 12) 'Sbox0Mm/Kg
        Workbooks(File).Worksheets(Sal).Cells(M, 15) =
Workbooks(File).Worksheets(BTL).Cells(j, 13) 'OxsolMm/Kg
    End For
End If

```

```

        Workbooks(File).Worksheets(Sal).Cells(M, 16) =
Workbooks(File).Worksheets(BTL).Cells(j, 14) 'PrDM
        Workbooks(File).Worksheets(Sal).Cells(M, 17) =
Workbooks(File).Worksheets(BTL).Cells(j, 15) 'T090C
        Workbooks(File).Worksheets(Sal).Cells(M, 18) =
Workbooks(File).Worksheets(BTL).Cells(j, 16) 'T190C
        Workbooks(File).Worksheets(Sal).Cells(M, 19) =
Workbooks(File).Worksheets(BTL).Cells(j, 17) 'C0S/m
        Workbooks(File).Worksheets(Sal).Cells(M, 20) =
Workbooks(File).Worksheets(BTL).Cells(j, 18) 'C1S/m
        Workbooks(File).Worksheets(Sal).Cells(M, 21) =
Workbooks(File).Worksheets(BTL).Cells(j, 19) 'Sbeox0V
        M = M + 1
    End If

```

```

    Next j
End If
Next i
End Sub

```

Function SolO2mLL(Temp, Salinity)

Rem : calculate oxygen % saturation @ theta

Rem H.E. GARCIA and L.I. GORDON, limnol. Oceanogr., 37(6), 1992,
1307-1312, Benson and Krause constants (ml/l)

```

    aox0 = 2.00907
    aox1 = 3.22014
    aox2 = 4.0501
    aox3 = 4.94457
    aox4 = -0.256847
    aox5 = 3.88767
    box0 = -0.00624523
    box1 = -0.00737614
    box2 = -0.010341
    box3 = -0.00817083
    cox0 = -4.88682E-07
    Ts = Log((298.15 - Temp) / (273.15 + Temp))
    Aox = aox0 + aox1 * Ts + aox2 * Ts ^ 2 + aox3 * Ts ^ 3 + aox4 * Ts ^ 4 + aox5
* Ts ^ 5
    Box = (box0 + box1 * Ts + box2 * Ts ^ 2 + box3 * Ts ^ 3) * Salinity
    Cox = Aox + Box + cox0 * Salinity ^ 2
    C4 = Exp(Cox) 'Solubility in ml/l
    SolO2mLL = C4
    'C4 = C4 * 1000 / 22.3916 ' O2 solubility in umol/l
End Function

```

Function SolO2μmolkg(Temp, Salinity)

'function [conc_O2] = O2sol(S,T)

'% O2sol Solubility of O2 in sea water

'%=====

=====
'% O2sol Version 1.1 4/4/2005

'% Author: Roberta C. Hamme (Scripps Inst of Oceanography)

'%

'% USAGE: concO2 = O2sol(S,T)

'%

'% DESCRIPTION:

'% Solubility (saturation) of oxygen (O2) in sea water

'% at 1-atm pressure of air including saturated water vapor

'%

'% INPUT: (if S and T are not singular they must have same dimensions)

'% S = salinity [PSS]

'% T = temperature [degree C]

'%

'% OUTPUT:

'% concO2 = solubility of O2 [umol/kg]

'%

'% AUTHOR: Roberta Hamme (rhamme@ucsd.edu)

'%

'% REFERENCE:

'% Hernan E. Garcia and Louis I. Gordon, 1992.

'% "Oxygen solubility in seawater: Better fitting equations"

'% Limnology and Oceanography, 37, pp. 1307-1312.

'%

'% DISCLAIMER:

'% This software is provided "as is" without warranty of any kind.

'%=====

=====
'% CALLER: general purpose

'% CALLEE: sw_dens_0.m

'%-----

'% Check input parameters

'%-----

'if nargin ~=2

' error('O2sol.m: Must pass 2 parameters')

'end %if

'% Check S,T dimensions and verify consistent

'[ms,ns] = size(S);

'[mt,nt] = size(T);

```

'
'% Check that T&S have the same shape or are singular
'if ((ms~=mt) | (ns~=nt)) & (ms+ns>2) & (mt+nt>2)
' error('O2sol: S & T must have same dimensions or be singular')
'end %if

'%-----
'% BEGIN
'%-----

'% convert T to scaled temperature
Temp_S = Log((298.15 - Temp) / (273.15 + Temp))

'% constants from Table 1 of Garcia & Gordon for the fit to Benson and Krause
(1984)
A0_O2 = 5.80871
A1_O2 = 3.20291
A2_O2 = 4.17887
A3_O2 = 5.10006
A4_O2 = -0.0986643
A5_O2 = 3.80369
B0_O2 = -0.00701577
B1_O2 = -0.00770028
B2_O2 = -0.0113864
B3_O2 = -0.00951519
C0_O2 = -2.75915E-07

' Corrected Eqn (8) of Garcia and Gordon 1992
Aox = A0_O2 + A1_O2 * Temp_S + A2_O2 * Temp_S ^ 2 + A3_O2 * Temp_S ^ 3 +
A4_O2 * Temp_S ^ 4 + A5_O2 * Temp_S ^ 5
Box = Salinity * (B0_O2 + B1_O2 * Temp_S + B2_O2 * Temp_S ^ 2 + B3_O2 *
Temp_S ^ 3)
Cox = Aox + Box + C0_O2 * Salinity ^ 2
SolO2μmolkg = Exp(Cox)

End Function

Function Density(Temp, Salinity)

'https://unesdoc.unesco.org/ark:/48223/pf0000188170
'Definitions cstes
'Density Constants
DSB0 = 0.824493
DSB1 = -0.0040899
DSB2 = 7.6438E-05

```

DSB3 = -8.2467E-07
DSB4 = 5.3875E-09

DSC0 = -0.00572466
DSC1 = 0.00010227
DSC2 = -1.6546E-06

DSD0 = 0.00048314

DSA0 = 999.842594
DSA1 = 0.06793952
DSA2 = -0.00909529
DSA3 = 0.0001001685
DSA4 = -1.120083E-06
DSA5 = 6.536332E-09

'Secant bulk modulus (K) of seawater : const

'Validity of EOS80 is valid for S=0 to 42; t=-2 to 40 °C, p=0 to 10000 db

DensityWaterTemp = DSA0 + DSA1 * Temp + DSA2 * Temp ^ 2 + DSA3 * Temp ^ 3 +
DSA4 * Temp ^ 4 + DSA5 * Temp ^ 5

DensityTemp = DensityWaterTemp + (DSB0 + DSB1 * Temp + DSB2 * Temp ^ 2 +
DSB3 * Temp ^ 3 + DSB4 * Temp ^ 4) * Salinity + (DSC0 + DSC1 * Temp + DSC2 * Temp ^
2) * Salinity ^ 1.5 + DSD0 * Salinity ^ 2

Density = DensityTemp / 1000

End Function

Function ThetaPot(Pressure, TIS, Salinity)

'Add references

'Definitions cste

'Density Constants

DSB0 = 0.824493
DSB1 = -0.0040899
DSB2 = 7.6438E-05
DSB3 = -8.2467E-07
DSB4 = 5.3875E-09

DSC0 = -0.00572466
DSC1 = 0.00010227
DSC2 = -1.6546E-06

DSD0 = 0.00048314

DSA0 = 999.842594
DSA1 = 0.06793952
DSA2 = -0.00909529
DSA3 = 0.0001001685
DSA4 = -1.120083E-06
DSA5 = 6.536332E-09

'Secant bulk modulus (K) of seawater : const

'Validity of EOS80 is valid for S=0 to 42; t=-2 to 40 °C, p=0 to 10000 db

KF0 = 54.6746

KF1 = -0.603459

KF2 = 0.011

KF3 = -6.17E-05

KG0 = 0.0794

KG1 = 0.0165

KG2 = -0.00053

KI0 = 0.00228

KI1 = -1.1E-05

KI2 = -1.61E-06

KJ0 = 0.000191

KM0 = -9.93E-07

KM1 = 2.08E-08

KM2 = 9.17E-10

KE0 = 19652.21

KE1 = 148.4206

KE2 = -2.327105

KE3 = 0.0136

KE4 = -5.16E-05

KH0 = 3.239908

KH1 = 0.00144

KH2 = 0.000116

KH3 = -5.78E-07

KK0 = 8.51E-05

KK1 = -6.12E-06

KK2 = 5.28E-08

'Adiabatic correction for theta pot

ATG0 = -2.1687E-16

ATG1 = 1.8676E-14

ATG2 = -4.6206E-13

ATG3 = 2.7759E-12

ATG4 = -1.1351E-10

ATG5 = -5.4481E-14

ATG6 = 8.733E-12

ATG7 = -6.7795E-10

ATG8 = 1.8741E-08

ATG9 = -4.2393E-08

ATG10 = 1.8932E-06

ATG11 = 6.6228E-10
ATG12 = -6.836E-08
ATG13 = 8.5258E-06
ATG14 = 3.5803E-05

Sal0 = 35
Pr0 = 0

TempIS = TIS
If Pressure = "" Then GoTo nexti
If TempIS = "" Then GoTo nexti
If Salinity = "" Then GoTo nexti
'Density pure water

DensityWaterIS = DSA0 + DSA1 * TempIS + DSA2 * TempIS ^ 2 + DSA3 * TempIS ^ 3 + DSA4 * TempIS ^ 4 + DSA5 * TempIS ^ 5

DensityIS = DensityWaterIS + (DSB0 + DSB1 * TempIS + DSB2 * TempIS ^ 2 + DSB3 * TempIS ^ 3 + DSB4 * TempIS ^ 4) * Salinity + (DSC0 + DSC1 * TempIS + DSC2 * TempIS ^ 2) * Salinity ^ 1.5 + DSD0 * Salinity ^ 2

'Calculation Theta

H = Pr0 - Pressure

ATG0_SP0T0 = (((ATG0 * TempIS + ATG1) * TempIS + ATG2) * Pressure + ((ATG3 * TempIS + ATG4) * (Salinity - Sal0) + ((ATG5 * TempIS + ATG6) * TempIS + ATG7) * TempIS + ATG8)) * Pressure + (ATG9 * TempIS + ATG10) * (Salinity - Sal0) + ((ATG11 * TempIS + ATG12) * TempIS + ATG13) * TempIS + ATG14

XK0 = H * ATG0_SP0T0

T1 = TempIS + 0.5 * XK0

Q0 = XK0

P1 = Pressure + 0.5 * H

AGT1_SP1T1 = (((ATG0 * T1 + ATG1) * T1 + ATG2) * P1 + ((ATG3 * T1 + ATG4) * (Salinity - Sal0) + ((ATG5 * T1 + ATG6) * T1 + ATG7) * T1 + ATG8)) * P1 + (ATG9 * T1 + ATG10) * (Salinity - Sal0) + ((ATG11 * T1 + ATG12) * T1 + ATG13) * T1 + ATG14

XK1 = H * AGT1_SP1T1

T2 = T1 + 0.29289322 * (XK1 - Q0)

Q1 = 0.58578644 * XK1 + 0.121320344 * Q0

ATG2_SP1T2 = (((ATG0 * T2 + ATG1) * T2 + ATG2) * P1 + ((ATG3 * T2 + ATG4) * (Salinity - Sal0) + ((ATG5 * T2 + ATG6) * T2 + ATG7) * T2 + ATG8)) * P1 + (ATG9 * T2 + ATG10) * (Salinity - Sal0) + ((ATG11 * T2 + ATG12) * T2 + ATG13) * T2 + ATG14

XK2 = H * ATG2_SP1T2

P2 = P1 + 0.5 * H

T3 = T2 + 1.707106781 * (XK2 - Q1)

ATG3_SP2T3 = (((ATG0 * T3 + ATG1) * T3 + ATG2) * P2 + ((ATG3 * T3 + ATG4) * (Salinity - Sal0) + ((ATG5 * T3 + ATG6) * T3 + ATG7) * T3 + ATG8)) * P2 + (ATG9 * T3 + ATG10) * (Salinity - Sal0) + ((ATG11 * T3 + ATG12) * T3 + ATG13) * T3 + ATG14

Q2 = 3.414213562 * XK2 - 4.121320344 * Q1

XK3 = ATG3_SP2T3 * H

Theta = (XK3 - 2 * Q2) / 6 + T3

ThetaPot = Theta

nexti:

End Function

Function Conductivity_Salinity(TIS, P, Conductivity)

'Conversion conductivity Salinity

'Cste definition

A1 = 2.07E-05

A2 = -6.37E-10

A3 = 3.989E-15

B1 = 0.03426

B2 = 0.0004464

B3 = 0.4215

B4 = -0.003107

C0 = 0.6766097

C1 = 0.0200564

C2 = 0.0001104259

C3 = -6.9698E-07

C4 = 1.0031E-09

Dim aa(6)

aa(0) = 0.008

aa(1) = -0.1692

aa(2) = 25.3851

aa(3) = 14.0941

aa(4) = -7.0261

aa(5) = 2.7081

Dim bb(6)

bb(0) = 0.0005

bb(1) = -0.0056

bb(2) = -0.0066

bb(3) = -0.0375

bb(4) = 0.0636

bb(5) = -0.0144

'P in decibar,T in deg C,C en S/m

'// C = conductivity S/m, T = temperature deg C ITPS-68, P = pressure in

decibars

If Conductivity <= 0# Then

Salinity = 0#

Else

Conductivity = 10# * Conductivity /* convert Siemens/meter to mmhos/cm */

R = Conductivity / 42.914

Value = 1 + B1 * TIS + B2 * TIS * TIS + B3 * R + B4 * R * TIS

RP = 1 + (P * (A1 + P * (A2 + P * A3))) / Value

Value = RP * (C0 + (TIS * (C1 + TIS * (C2 + TIS * (C3 + TIS * C4))))))

RT = R / Value

If (RT <= 0#) Then RT = 1E-06

sum1 = 0

sum2 = 0

```

        For ii = 0 To 5
            Temp = RT ^ (ii / 2#)
            sum1 = sum1 + aa(ii) * Temp
            sum2 = sum2 + bb(ii) * Temp
        Next ii
        Value = 1# + 0.0162 * (TIS - 15#)
        Salinity = sum1 + sum2 * (TIS - 15#) / Value
        Conductivity_Salinity = Salinity
    End If
End Function

```

```

Sub FILTER_OUT()

```

```

File = "IAPSO_BTL_Merged.xlsm"
SBE43_O2W = "O2WSBE43_Ox0_v02"
OUT = "Removed"
Start_i = 2
End_i = 1136
Start_k = 4
End_k = 1136
P = 1
MyLoop = 1
Do_It_Again:
    For i = Start_i To End_i
        flag = Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 26)
        If flag = "OUT" Then
            'Select Line
            MyRange = "(d" & i & ":Y" & i & ")"
            MyFormula = MyRange
            Range(MyFormula).Select
            'Erase line
            Selection.Cut
            Worksheets(Removed).Activate'
            Workbooks(File).Worksheets(OUT).Cells(P, 1) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 1)
            '
            Workbooks(File).Worksheets(OUT).Cells(P, 2) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 2)
            '
            Workbooks(File).Worksheets(OUT).Cells(P, 1) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 3)

            Selection.ClearContents
        End If
    Next i
    SolverOk SetCell:="$V$2", MaxMinVal:=2, ValueOf:=0, ByChange:="$AA$1:$AA$3",
Engine:=3, EngineDesc:="Evolutionary"

```

```

'SolverSolve
SolverSolve UserFinish:=True
SolverFinish KeepFinal:=1
For j = Start_i To End_i
    Flag2 = Workbooks(File).Worksheets(SBE43_O2W).Cells(j, 26)
    If Flag2 = "OUT" Then
        MyLoop = MyLoop + 1
        GoTo Do_It_Again
    End If
Next j
Workbooks(File).Worksheets(SBE43_O2W).Cells(32, 28) = MyLoop
For k = Start_k To End_k
    If Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 4) = "" Then
        Workbooks(File).Worksheets(OUT).Cells(P, 1) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 1)
        Workbooks(File).Worksheets(OUT).Cells(P, 2) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 2)
        Workbooks(File).Worksheets(OUT).Cells(P, 3) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 3)
        P = P + 1
    End If
Next k
End Sub

```

```

Sub Salinity_FILTER_OUT()

```

```

    'Salinity_SBE04_00_1
    ,

```

```

File = "IAPSO_BTL_Merged.xlsm"
SBE43_O2W = "Salinity_SBE04_00_1"
OUT = "Removed_Salinity"
Start_i = 4
End_i = 1151
Start_k = 4
End_k = 1151
P = 1
MyLoop = 1
Flag2 = ""
Do_It_Again:
    For i = Start_i To End_i
        flag = Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 13)
        If flag = "OUT" Then
            'Select Line
            MyRange = "(d" & i & ":M" & i & ")"

```

```

        MyFormula = MyRange
        Range(MyFormula).Select
        'Erase line
    '
        Selection.Cut
    '
        Worksheets(Removed).Activate'
    '
        Workbooks(File).Worksheets(OUT).Cells(P, 1) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 1)
    '
        Workbooks(File).Worksheets(OUT).Cells(P, 2) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 2)
    '
        Workbooks(File).Worksheets(OUT).Cells(P, 1) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(i, 3)

        Selection.ClearContents
    End If
Next i
    SolverOk SetCell:="$V$2", MaxMinVal:=2, ValueOf:=0,
ByChange:="$NN1$1:$NN$3", Engine:=3, EngineDesc:="Evolutionary"
    'SolverSolve
    SolverSolve UserFinish:=True
    SolverFinish KeepFinal:=1
    For j = Start_i To End_i
        Flag2 = Workbooks(File).Worksheets(SBE43_O2W).Cells(j, 26)
        If Flag2 = "OUT" Then
            MyLoop = MyLoop + 1
            GoTo Do_It_Again
        End If
    Next j
    Workbooks(File).Worksheets(SBE43_O2W).Cells(20, 15) = MyLoop
    For k = Start_k To End_k
        If Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 4) = "" Then
            Workbooks(File).Worksheets(OUT).Cells(P, 1) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 1)
            Workbooks(File).Worksheets(OUT).Cells(P, 2) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 2)
            Workbooks(File).Worksheets(OUT).Cells(P, 3) =
Workbooks(File).Worksheets(SBE43_O2W).Cells(k, 3)
            P = P + 1
        End If
    Next k
End Sub

```